

WISH IN IVTLatex

A Quantitative ImmunoDiagnostic Method

Product Features & Customize Service

The turbidimetric assay is based on the agglutination reaction between latex particles coated with antibody and the antigen in solution.

We support your application needs by offering custom specifications to meet your precise requirements.

ITEMS

1. Anti-CCP (anti-cyclic citrullinated peptide antibody)
2. PCIII (Precollagen Type III)
3. CIV (Collagen Type IV)
4. HPAb (Helicobacter Pylori Antibody)

Coming Soon

1. IL-6 (Interleukin-6)
2. G17 (Gastrin-17)
3. MMP3 (Matrix Metalloproteinase 3)
4. KL-6 (Krebs von den Lungen 6)